




2606KIME.332L
KC-AD8-01 ACD856
DMSO, H-1
NUMEGA 10-19-2026




2606KIME.332L
KC-AD8-01 ACD856
DMSO, H-1
NUMEGA 10-19-2026




2606KIME.332L
KC-AD8-01 ACD856
DMSO, H-1
NUMEGA 10-19-2026




2606KIME.332T
KC-AD8-01 ACD856
DMSO, C-13
NUMEGA 10-19-2026




2606KIME.332T
KC-AD8-01 ACD856
DMSO, C-13
NUMEGA 10-19-2026



2606KIME.332T
KC-AD8-01 ACD856
DMSO, C-13
NUMEGA 10-19-2026



2606KIME.332D
KC-AD8-01 ACD856
DMSO, QH-1
NUMEGA 10-15-2026



October 20, 2026

NuMega Reference #
2606KIME.332D

Company: Kimera Chem

Quantitative Analysis Report

Bruker 500 MHz, AVII-500 and AV-500

Sample ID	KC-AD8-01 ACD856
Standard	Ethylene Carbonate
Weight of Standard (wt _{st})	1.976
Weight of Sample (wt _{sa})	12.972
Assumed Sample Molecular Weight (MW _{sa})	387.4
Standard Molecular Weight (MW _{st})	88.06
Integration of Sample (I _{sa})	36.55
Integration of Standard (I _{st})	100
# of Equivalent Nuclei in Standard (n _{st})	4
# of Equivalent Nuclei in Sample (n _{sa})	1
Assigned Purity of Standard % (P)	99.71
Purity %	97.7

Formula Used for Calculation:

$$\text{Purity \%} = \frac{wt_{st}}{wt_{sa}} \times \frac{MW_{sa}}{MW_{st}} \times \frac{I_{sa}}{I_{st}} \times \frac{n_{st}}{n_{sa}} \times P$$